

PrograBase 2

Control

Manage

Protocol


Download

The multifunctional software for your Programat®



PrograBase 2

The multifunctional software for your Programat®



CONTROL

Control your Programat ceramic furnace directly from your computer. Relevant information, such as program number, firing status and current firing temperature are indicated on the screen.

MANAGE

You can easily manage your firing programs via your computer and subsequently save the data on the furnace

and/or your computer. You can change the firing programs as you wish, create new programs and copy existing data from one furnace to another.

PROTOCOL

Protocol the ceramic firing cycles you have conducted by means of PrograBase 2. The protocolled firing cycles can be supplemented by additional data, such as patient information, order number or material information. You can save and print the protocols if you like.



DOWNLOAD

Our homepage www.ivoclar-vivadent.com can be accessed with one mouse click with the PrograBase 2 software to download the latest equipment software, Operating Instructions and firing tables for free.

System requirements

- Operating system: Windows XP® (including Microsoft .NET Framework 2.0*), Windows Vista® or Windows 7
- Screen resolution 1024 x 768
- Colours: min. 256
- RAM: min. 256MB
- CPU: 1GHz
- Min. 11 MB free storage space on hard drive
- Internet access

*) If Microsoft .NET Framework 2.0 is not available on your computer, it will be installed together with the software. Required free storage space about 300 MB.

Supported furnaces

The PrograBase 2 supports the following furnaces:

- Programat P300/G2
- Programat P500/G2 + P500
- Programat P700/G2 + P700
- Programat EP 3000/G2 + EP 3000
- Programat EP 5000/G2 + EP 5000
- Programat CS
- Programat S1

Required accessories

USB download cable (included in the delivery form of the furnace)



Programat P300/G2



Programat P500/G2



Programat P700/G2



Programat EP 3000/G2



Programat EP 5000/G2



Programat CS



Programat S1

Descriptions and data constitute no warranty of attributes.

Printed in Austria
© Ivoclar Vivadent AG
614929/0311/e

Ivoclar Vivadent AG
Bendererstr. 2
FL-9494 Schaan
Principality of Liechtenstein
Phone +423 / 235 35 35
Fax +423 / 235 33 60
www.ivoclarvivadent.com

ivoclar
vivadent
passion vision innovation